

RINSHO HINYOKIKA

(JAPANESE JOURNAL OF CLINICAL UROLOGY)

Advisory Board : **Toyohei Machida MD, Akimi Ogawa MD, Masao Akimoto MD, Masaru Murai MD**

Editorial Board : **Kenjiro Kohri MD, Tomoaki Fujioka MD, Mototsugu Oya MD**

CONTENTS

Vol 66 No 6 May 2012

SPECIAL THEME : Selection and directions of the energy devices

How to use monopolar electro surgical device	Yukio Naya	365
Repertoire of electrosurgical devices with bipolar mode in urologic surgery	Kazunobu Shinoda, et al	371
Properties and surgical usage of a microwave tissue coagulator	Satoshi Anai, et al	379
Ultrasonically activated device	Hidefumi Kinoshita, et al	385
How to use vessel sealing system	Daisuke Nagata	391
Argon beam coagulator	Nobuo Tsuru, et al	399
How to choose and use radiofrequency incision/coagulation device	So Omori	405
Choice and usage of Ho : YAG laser	Kogenta Nakamura	411
Utilization of mechanical stapler in urologic surgery	Ken Takeda, et al	417

CASE REPORTS

Renal pelvic carcinoma producing granulocyte colony-stimulating factor (gcsf) : a case report	Naoya Fujikawa, et al	429
---	-----------------------	-----